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
Original operating instructions

First steps

What do these operating instructions apply to?

These operating instructions apply to the following devices, which differ in terms of their range of functions or equipment features:

FN 6864/23 MediFlex for room radio
 FN 6064/23 MediFlex for room bus

Identifying the device

Make sure that these operating instructions belong to your device. These operating instructions do not apply to other devices, even if they look similar or appear to be of the same design.

Unique feature: Type plate ("x" = any number)

Position of the type plate:


On the electronics unit

Checking the scope of delivery

Check that the scope of delivery is complete and in perfect condition. In case of missing parts or damage: Do not connect the device and do not put it into operation! Make a complaint to the supplier.

The scope of delivery consists of:

- 1 Paddle with proximity sensor and electronics unit
- 2 Articulated arm
- 3 CR2450 battery (installed ready for use)
- 4 Download the latest version of the operating instructions at <https://www.tetronik-kommunikationstechnik.de/downloads/original-betriebsanleitungen/>. Search term: FN number of your device.

What to do with the packaging?

If the device or parts of it are to be stored before or after use: Keep the packaging.

If nothing is to be stored:

Dispose of the packaging in accordance with regional regulations. Do not dispose of the packaging in household waste!

Using the operating instructions

Read first!

Read these operating instructions in full. Keep the operating instructions in a safe place. If you have any questions, please contact the manufacturer.

Interpreting signal words

Pay attention to signal words that warn of dangers. Follow the instructions for risk avoidance.

Signal word	Meaning
DANGER	Warns of imminent danger to life or danger of serious injury.
WARNING	Warns of possible imminent danger to life or danger of serious injury.
CAUTION	Warns of possible imminent danger of minor or slight injuries.
SAFETY INSTRUCTIONS	Warns of possible material damage.
NOTICE	Indicates additional information.


Information about the device

Intended use

MediFlex enables residents or patients with severe motor impairments to trigger a call by bringing the palm of the hand close to the sensor or by slightly bending the paddle.

Any other use is contrary to the intended use and may lead to malfunctions and damage to this and connected devices. Disregarding these operating instructions also constitutes improper use.

Arrangement in the call system


Features and properties

- Call triggering by approaching or bending
- Wide range of fastening options using clamping jaws
- Individual alignment via articulated arm
- Acoustic and visual call display


Controls and displays

The device has the following controls and displays:


On the paddle:

- 1 Foldable surface
- 2 Call display
- 3 Proximity sensor

On the electronics unit:

- 4 Battery compartment
- 5 Holder for the articulated arm
- 6 DIP switch for configuration

On the articulated arm:

- 7 Paddle attachment
- 8 Fixing toggle
- 9 Clamping jaws

Technical data

Physical data

Paddle dimensions (W x H x D)	320 x 90 x 28 mm
Length of articulated arm (extended)	max. 270 mm
Weight incl. articulated arm	approx. 490 g
Paddle colour	Orange

Electrical data

Supply voltage	Battery CR2450
Battery life	Approx. 3 years
Operating frequency	433.92 MHz
Protection class	IP 44
Transmitter range	approx. 10 metres

Ambient conditions

Permissible ambient temperature	0...+55°C
Permissible storage temperature	-10...+60°C
Relative humidity	10...95%

The device is only suitable for use in dry indoor areas. It must not be exposed to moisture.

The device is only suitable for use in shielded rooms (e.g. X-ray rooms) if the associated receiver is located in the same room.

The device must not be used in potentially explosive atmospheres.


CE The declaration of conformity can be requested from the manufacturer. See last page for contact details.

Normative information

- Planning and testing of the call system and certification of its functionality must be carried out by a "specialist planner for call systems".
- This device is not a medical device within the meaning of Directive 93/42/EEC.
- The electrical safety of the call system is guaranteed by system separation.

Instructions for technical staff

Procuring the required accessories

- 1 FN 6864/00
IMA receiver for room radio
- 2 Receiver for wireless handheld button
FN 6064/10 7-pole DIN
FN 6064/50 8-pole mini-DIN

Putting the device into operation

Ensure that the following requirements are met:

- The operating instructions for the associated receiver are to hand.
- The associated receiver operates on the same frequency as the call pad.
- The associated receiver is ready for use and is within reach of the call pad.


Configuring the device

The required configuration of the device depends on the condition of the resident or patient. Coordinate with the staff in charge which functions are required or desired.

The device is switched off on delivery.

- 1 Switch on the device.
- 2 Should the device beep 5 times when a call is triggered?


ON Switch on the loudspeaker.

OFF Switch off the loudspeaker.

- 3 Should the device react when approached by a hand? *)


ON Switch on the proximity sensor.

OFF Switch off the proximity sensor.

- 4 Should the device react to bending? *)


ON Switch on the bend detection.

OFF Switch off the bend detection.


*) Functions 3 and 4 can be activated individually or together as required.

Registering the device on the receiver

- 1 Trigger a call.
 - The proximity sensor flashes 5 times in red (call indicator).
 - When the loudspeaker is switched on: The device beeps 5 times.
- 2 Continue the registration process as described in the operating instructions for the corresponding receiver. MediFlex is then ready for use.

Replacing the battery

Battery type: CR2450

The battery is located in the battery compartment of the electronics unit.

- 1 Open the battery compartment.
- 2 Lever out the battery.
- 3 Insert the new battery.
Ensure correct polarity (positive pole at the top).
- 4 Close the battery compartment.

De-registering the device

If the device is no longer needed on the associated receiver:

- 1 Switch off the device.
- 2 Continue the de-registration process as described in the operating instructions for the corresponding receiver.
MediFlex is then de-registered.

Servicing and maintaining the device

Ensure that this device and the entire call system are fail-safe at all times. Follow the regulation DIN VDE 0834-1:2016-06 (11.2 to 11.4).

Repairing the device

The device fulfils safety-relevant tasks.

There is a risk of malfunction if the device is repaired incorrectly.


Only have repairs carried out by the manufacturer.

Instructions for operating staff


Attaching the device

The device is mounted on a bed or wheelchair, for example. Installation is possible wherever there is a suitable rod for the articulated arm.


The articulated arm must be attached in such a way that the paddle can be positioned conveniently for the resident or patient. The user must be able to reach the device effortlessly in accordance with their mobility.

- 1 Loosen the fixing toggle. This allows all joints and rotating parts to move.
- 2 Attach the paddle to the articulated arm.
- 3 Place the clamping jaws on a rod.
- 4 Close the clamping jaws until they firmly grip the rod.


- 5 The articulated arm can be moved freely in all directions. Find the desired position for the paddle and hold it there.

- 6 Tighten the fixing toggle firmly. This locks all joints and rotating parts simultaneously and the articulated arm is fixed.


Triggering a call

To trigger a call, proceed as follows depending on the configuration of the device:

Cover the proximity sensor with the palm of your hand.

Or:

Bend the paddle.

The call indicator of the proximity sensor flashes 5 times (reassurance light).

When the loudspeaker is switched on: The device beeps 5 times.

Teaching residents or patients

Explain all the functions that are important to the resident or patient:

- Triggering a call

Instructions for auxiliary staff

Cleaning

Clean the surfaces with a damp, lint-free cloth. Do not use any aggressive or abrasive cleaning agents.

Disinfecting

Disinfect the surface of the device with a commercially available surface disinfectant according to the instructions.

SAFETY INSTRUCTIONS

The device contains heat-sensitive components.
There is a risk of destruction when exposed to heat.
Do not sterilise the device thermally.

Procuring consumables

SAFETY INSTRUCTIONS

The device fulfils safety-relevant tasks.
There is a risk of malfunction if foreign material is used.
Only use original material from the manufacturer.

The following consumables are available:

Battery CR2450

Storing the device when not in use

Store the device in its original packaging in a suitable storage location. For the requirements for the storage location, see "Technical data".

Disposing of the device

At the end of the service life:

Observe the regional disposal regulations. Have the device recycled properly. Do not dispose of the device with household waste!